

Catalog Number

9502S-E-0AV

Vials, Screw Top, Glass. Amber, 2ml, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Borosilicate Glass, amberized. Features a wide mouth for easy filling. Lowest cost in class.

Amber vials are often used for light sensitive compounds and considered low actinic.

Price
\$33.19

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$33.19