

Catalog Number

9502S-E-10C-B

Caps, Screw Top, with Soft-Guard, White Silicone Rubber / Red PTFE Septa, bonded in "knurled", polypropylene, Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. Soft-Guard Brand. 100/PK



Additional Info

This is our most popular screw cap for autosampler vials.

Soft-Guard™ Septa are easier to penetrate by autosampler needles, and re-seals better than other Caps on the market. This provides longer instrument up-time, fewer missed injections and more consistent data. Use for HPLC, GC, LCMS or GCMS. Click [HERE](#) for more information about Soft-Guard™ Septa Technology.

The Durometer of the Septa is 40-50. The thickness of the Septa is 0.040" (+/- .003")

Click [HERE](#) for case quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$36.79

Printed On 2024/05/01