

## Catalog Number

9502S-E-10C-N

**Caps, Screw Top**, with Soft-Guard, White Silicone Rubber / Red PTFE Septa, bonded in "knurled", polypropylene, Natural Translucent Caps. For use with 9-425mm thread, 12x32mm autosampler vials. Soft-Guard



### Additional Info

Soft-Guard™ Septa are easier to penetrate by autosampler needles, and re-seals better than other Caps on the market. This provides longer instrument up-time, fewer missed injections and more consistent data. Use for HPLC, GC, LCMS or GCMS. Click [HERE](#) for more information about Soft-Guard™ Septa Technology.

Click [HERE](#) for case quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$36.79**

Printed On 2024/05/02