

Catalog Number**9502S-E-10C-P-M**

Caps, Screw Top, with Soft-Guard, White Silicone Rubber / Red PTFE Septa, bonded in "knurled", polypropylene, Purple Caps. For use with 9-425mm thread, 12x32mm autosampler vials. Soft-Guard Brand. 1000/CS

**Additional Info**

Soft-Guard™ Septa are easier to penetrate by autosampler needles, and re-seals better than other Caps on the market. This provides longer instrument up-time, fewer missed injections and more consistent data. Use for HPLC, GC, LCMS or GCMS.

Click [HERE](#) for more information about Soft-Guard™ Septa Technology.

Click [HERE](#) for smaller pack quantities.

Price**\$357.79**

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$357.79