

Catalog Number

9502S-E-15K-M

Caps, Screw Top, with White Silicone Rubber / Natural PTFE Septa, fitted in "knurled", polypropylene, Black Caps. For use with 9-425mm thread, 12x32mm autosampler vials. BASIK Brand. 1000/CS.



Additional Info

This is our most economical, lowest cost Screw Cap for Autosampler Vials.

Natural PTFE is as equally inert and durable as Red PTFE, but no dyes have been added. Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$274.16

Printed On 2024/02/25