


## Catalog Number

9502S-E-15Y

**Caps, Screw Top**, with White Silicone Rubber / Natural PTFE Septa, fitted in "knurled", polypropylene, Yellow Caps. For use with 9-425mm thread, 12x32mm autosampler vials. BASIK Brand. 100/PK.



### Additional Info

This is our most economical, lowest cost Screw Cap for Autosampler Vials.

Natural PTFE is as equally inert and durable as Red PTFE, but no dyes have been added. Click [HERE](#) for case quantities.

**Price**  
**\$30.64**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$30.64**