

Catalog Number

9502S-E-34B

Caps, Screw Top, with Soft-Guard, White Silicone Rubber / Red PTFE Pre-Slit (Y Shape) Septa, fitted in "knurled", polypropylene, Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. Soft-Guard



Additional Info

Pre-slit for easy penetration by autosampler needles. Recommended for Waters® Alliance & Acquity instruments. Y Shape Pre-Slit provides a greater target area for needles and are easier to penetrate but do not reseal as well as single slit septa.

Soft-Guard™ Septa are easier to penetrate by autosampler needles, and re-seals better than other Caps on the market. This provides longer instrument up-time, fewer missed injections and more consistent data. Use for HPLC, GC, LCMS or GCMS. Click [HERE](#) for more information about Soft-Guard™ Septa Technology.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$46.71

Printed On 2024/05/06