


Catalog Number

9502S-E-35R

Caps, Screw Top, with White Silicone Rubber / Natural PTFE Pre-Slit Septa, fitted in "knurled", polypropylene, Red Caps. For use with 9-425mm thread, 12x32mm autosampler vials. BASIK Brand. 100/PK.



Additional Info

This is our lowest cost Screw Cap with Pre-Slit Septa for Autosampler Vials.

Natural PTFE is equally inert and as durable as Red PTFE, but dyes have been added to it. Normally recommended for Waters® Alliance and Acquity instruments.

Click [HERE](#) for case quantities.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$34.68

Printed On 2024/02/25