

## Catalog Number

9502S-E-3DAP

**Vial & Cap Kit** Includes 100 EA 2ml, Screw Top, Amber Deactivated Autosampler Vials and 100 EA matching Screw Caps with White Silicone Rubber / Red PTFE Pre-Slit, Soft-Guard Septa, bonded in the Blue



### Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps. The kit packaging will open and re-close making it useful for keeping vial and cap types matched, organized and clean during use.

Click on links below for more Specifications for the 9mm Vials and Caps:

Supplied with Silanized / Deactivated Vials: Cat. No [9502S-E-0DAV](#) and Caps: Cat No. [9502S-E-30C-B](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$91.34**

Printed On 2024/05/05