

## Catalog Number

9502S-T-3CP

**Vial & Cap Kit** Includes 100 EA 2ml, Screw Top, Clear Autosampler Vials and 100 EA matching Screw Caps with Silicone Rubber / PTFE, ultra-pure Pre-Slit Septa, fitted in the Light Blue Caps. BASIK Brand Easy



### Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps. The kit packaging will open and re-close making it useful for keeping vial and cap types matched, organized and clean during use.

Click on links below for more Specifications for the 9mm Vials and Caps:

Supplied with Vials Cat. No [9502S-E-0CV](#) and Caps: Cat No. [9509S-30C-B](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$58.26**

Printed On 2024/04/28