

Catalog Number

9502SB-03

Storage Box for 2ml Autosampler Vials. Pink, chemical resistant polypropylene holds 81 vials. Complete with translucent plastic lid. 1 EA.



Additional Info

Excellent for standing position on the lab bench or during transport and ideal for space saving storage in places such as refrigerators.

- Lid prevents condensation on the caps minimizing contamination during storage.
- Usable temperature range -80C to +100C.
- Alpha numeric indexing coded on bottom of wells and along the sides of the box for clear sample identification.
- Box outer dimensions are 130mm wide by 130mm deep by 45mm high.
- Box and lid are Autoclavable.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$43.60

Printed On 2024/05/18