

Catalog Number

9503S-0DCV

Vials, Screw Top, Glass. Clear, 2ml, Deactivated, 10-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial.

MICROSOLV Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type N51A, Deactivated / Silanized pharmaceutical grade, Borosilicate Glass.

Features a wide mouth opening for easy filling and ND10 threads.

Price
\$59.13

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$59.13