

Catalog Number

9503S-10C-K-M

Caps, Screw Top, with Silicone Rubber / PTFE Septa, bonded in wide mouth, "knurled", polypropylene, Black Caps. For use with 10-425mm threads, 12x32mm autosampler vials. BASIK Brand. 1000/CS.



Additional Info

Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$522.26

Printed On 2024/02/29