

Catalog Number

9503S-30C-K

Caps, Screw Top, with Silicone Rubber / PTFE Pre-Slit Septa, bonded in wide mouth, "knurled", polypropylene, Black Caps. For use with 10-425mm threads, 12x32mm autosampler vials. BASIK Brand. 100/PK.



Additional Info

Pre-slit for easy penetration by autosampler needles.

Click [HERE](#) for case quantities.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$69.46

Printed On 2024/08/06

Additional Info

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$69.46

Printed On 2024/08/06