

## Catalog Number

9503S-30C-K-M

**Caps, Screw Top**, with Silicone Rubber / PTFE Pre-Slit Septa, bonded in wide mouth, "knurled", polypropylene, Black Caps. For use with 10-425mm threads, 12x32mm autosampler vials. BASIK Brand. 1000/CS.



### Additional Info

Pre-slit for easy penetration by autosampler needles.

Click [HERE](#) for smaller pack quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$624.43**

Printed On 2024/04/29