

Catalog Number

9503S-35C-K-M

Caps, Screw Top, with Silicone Rubber / PTFE Pre-Slit, 40mil Septa, bonded in wide mouth, "knurled", polypropylene, Black Caps. For use with 10-425mm threads, 12x32mm autosampler vials. BASIK Brand.



Additional Info

Pre-slit for easy penetration by autosampler needles.

This Septa is slightly thinner than most Septa to enhance the ease or needle penetration. 40mil equals 1.016mm.

Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$491.76

Printed On 2024/05/03