


Catalog Number

9503S-E-11C-LB

Caps, Screw Top, with 1.35mm thick, Silicone Rubber / PTFE Septa fitted in "knurled", polypropylene, Light Blue Caps. For use with 10-425mm thread, 12x32mm wide mouth autosampler vials. BASIK Brand.



MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$32.38

Printed On 2024/05/24