

## Catalog Number

9503S-E-WCV

**Vials, Screw Top, Glass.** Clear, 2ml, with Write-On Patch and fill lines, 10-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial. BASIK Brand. 100/PK.



### Additional Info

First hydrolytic, Class A, Borosilicate Glass. Features a wide mouth for easy filling. Low cost for budget minded laboratories. Write-On Patch and fill lines are epoxy paint for reduced smearing with solvents.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$33.25**

Printed On 2024/05/02