

## Catalog Number

## 9503S-GR

**Caps, Screw Top**, with Silicone Rubber / PTFE Septa fitted in "knurled", polypropylene, Black Caps. For use with 10-425mm thread, 12x32mm wide mouth autosampler vials. Compare to Shimadzu Caps. BASIK



### Additional Info

Flange at the bottom of the Cap allows for the vial and cap to be suspended in older Shimadzu autosamplers.

Click [HERE](#) for case quantities.

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$31.93**

Printed On 2024/02/25