

## Catalog Number

9503S-WCV

**Vials, Screw Top, Glass. Clear, 2ml, with Write-On Patch and fill lines, 10-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. MICROSOLV Brand. 100/PK.**



### Additional Info

First hydrolytic, Class A, Type N51A pharmaceutical grade, Borosilicate Glass. Features a wide mouth opening for easy filling and ND10 threads. Write-On Patch and fill lines are epoxy paint for reduced smearing.

**Price**  
**\$35.95**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$35.95**