


## Catalog Number

95040-SB1

**Storage Box** for 20ml EPA Style Vials. Violet, chemical resistant polypropylene holds 16 vials. Complete with translucent plastic lid. 1 EA.



### Additional Info

Excellent for standing position on the lab bench or during transport and ideal for space saving storage in places such as refrigerators.

- Lid prevents condensation on the caps minimizing contamination during storage.
- Usable temperature range -80C to +100C.
- Box outer dimensions are 130mm wide by 130mm deep by 102mm high.
- Box and lid are Autoclavable.

**Price**  
**\$68.95**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$68.95**