

Catalog Number

9504S-00C-BK-M

Caps, Screw Top, with an Open Hole. "knurled", polypropylene, Black Caps. Septa are sold separately and must be inserted manually. For use with 13-425mm thread, 15x45mm autosampler vials. MicroSolv Brand.



MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$134.30

Printed On 2024/04/29