

Catalog Number

9504S-01ND

Inserts, Low Volume for Vials, Glass. Clear, Deactivated, 250ul fill volume with a conical, precision point and 6x38mm outer dimensions. For use in 4ml autosampler vials. AQ Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type 33, Deactivated / Silanized Borosilicate Glass. Features a mandrel manufactured, precision point conical tip for maximum recovery of sample and internal consistency from insert to insert. Residual volume is 3ul with most autosamplers with a 150ul useable volume. This product is manufactured in an ISO 9001:2015 Certified facility.

Springs may be required and are sold separately. Use with 13mm Screw Top 15x45mm autosampler vials.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$226.74

Printed On 2024/05/24