

Catalog Number

9504S-0AV

Vials, Screw Top, Glass. Amber, 4ml, 13-425mm threads, 15x45mm
outer dimensions. For use as an autosampler vial or storage bottle.
MicroSolv Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type N51A pharmaceutical grade, Borosilicate Glass, amberized.
Features ND13 threads.

Amber vials are often used for light sensitive compounds and considered low actinic.

Use these vials with Waters WISP® and older Shimadzu® autosamplers.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$51.41

Printed On 2024/05/06