

## Catalog Number

## 9504S-30C-BK-M

**Caps, Screw Top, with Silicone Rubber / PTFE Pre-Slit Septa, bonded in "knurled", polypropylene, Black Caps. For use with 13-425mm thread, 15x45mm autosampler vials. MicroSolv Brand. 1000/CS.**



### Additional Info

Pre-slit for easy penetration by autosampler needles. Recommended for older Waters® instruments. Use with 15x45mm, 4ml vials.

Click [HERE](#) for smaller pack quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$681.86**

Printed On 2024/02/25