

Catalog Number

9504S-80C-BK

Caps, Screw Top, Solid, with a Red Orange Natural Rubber/Transparent TFE liner, in "knurled", polypropylene, Black Caps. For use with 13-425mm thread, 15x45mm autosampler vials. MicroSolv Brand.



Additional Info

This Solid Cap turns 15x45mm, 4ml autosampler vials into storage bottles. Rubber and TFE liner provide an inert seal.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$53.80

Printed On 2024/07/27