

Catalog Number

9504S-80C-BK

Caps, Screw Top, Solid, with a Red Orange Natural Rubber/Transparent TFE liner, in "knurled", polypropylene, Black Caps. For use with 13-425mm thread, 15x45mm autosampler vials. MicroSolv Brand.



Additional Info

This Solid Cap turns 15x45mm, 4ml autosampler vials into storage bottles. Rubber and TFE liner provide an inert seal.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$53.80

Printed On 2024/02/25