

## Catalog Number

9504S-E-05CP

Vial & Cap Kit Includes 100 EA 4ml, Screw Top, Clear Autosampler Vials and 100 EA matching Screw Caps with Silicone Rubber / PTFE Septa, fitted in the Black Caps, BASIK Brand Easy Purchase Pack. 1/PK.



### Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps.

Click on links below for more Specifications for the 13mm Thread Vials and Caps:

Supplied with Vials: Cat. No [9504S-E-0CV](#) and Caps: Cat No. [9504S-05C-BLK](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$60.76**

Printed On 2024/05/05