

Catalog Number

9504S-E-WCV

Vials, Screw Top, Glass. Clear, 4ml, with Write-On Patch and graduations, 13-425mm threads, 15x45mm outer dimensions. For use as an autosampler vial. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Borosilicate Glass. Features a wide mouth for easy filling. Low cost for budget minded laboratories. Write-On Patch and fill lines are epoxy paint for reduced smearing with solvents.

Price
\$40.54

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$40.54