

Catalog Number

9504S-WCV

Vials, Screw Top, Glass. Clear, 4ml, with Write-On Patch and fill lines, 13-425mm threads, 15x45mm outer dimensions. For use as an autosampler vial or storage bottle. MicroSolv Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type N51A pharmaceutical grade, Borosilicate Glass. Features ND13 threads. Write-On Patch and fill lines are epoxy paint for reduced smearing with solvents.

Use these vials with Waters WISP® and older Shimadzu® autosamplers.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$50.04

Printed On 2024/06/20