

Catalog Number

9504SB-03

Storage Box for 4ml Autosampler Vials. Red, chemical resistant polypropylene holds 49 vials. Complete with translucent plastic lid. 1 EA.



Additional Info

Excellent for standing position on the lab bench or during transport and ideal for space saving storage in places such as refrigerators.

- Lid prevents condensation on the caps minimizing contamination during storage.
- Usable temperature range -80C to +100C.
- Alpha numeric indexing coded on bottom of wells and along the sides of the box for clear sample identification.
- Box outer dimensions are 130mm wide by 130mm deep by 52mm high.

Price
\$63.82

Additional Info

- Box and lid are Autoclavable.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$63.82