

Catalog Number

9508C-0RB

Vials, Crimp Top, Glass. Clear, 0.3ml. An 8mm Serum Finish with 5.5x31.5mm outer dimensions and a round bottom. For use as an autosampler vial. MicroSolv Brand. (10 packs of 100/PK). 1000/CS.



Additional Info

First hydrolytic, Class A, pharmaceutical grade, Borosilicate Glass.

Use with 8mm aluminum crimp caps, 8mm or 9mm push caps. These round bottom vials have serum finish tops that provide the best seal for chromatography, when used properly with correct aluminum crimp caps.

These MicroVials have a round bottom outer shape and often need a holder / adapter to be used in an autosampler and will not stand by themselves.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$745.20

Printed On 2024/05/04