

## Catalog Number

9508C-30T-S

**Caps, Crimp Top**, with White Silicone Rubber / Red PTFE Pre-Slit Septa, Silver Caps, with a 4mm center hole. For use with 8mm Vials. MicroSolv Brand. 100/PK.



### Additional Info

These aluminum crimp caps are 8mm in diameter with a 4mm center hole and have a Pre-Slit, Ultra-Clean, fitted silicone rubber and PTFE septum that is 1.3mm thick. Use these caps with 8mm Serum Finish, Crimp Top Autosampler Vials.

These are Pre-Slit to make it easier for the instrument needle to penetrate the septum and offers reduced resistance. *It may not reseal as well as other Non-Slit Septa.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$82.89**

Printed On 2024/05/05