

Catalog Number

9508C-FB08-CV

Vials, Crimp Top, Glass. Clear, 0.8ml. An 8mm Serum Finish with 8.2x30mm outer dimensions and a flat bottom. For use as an autosampler vial. MicroSolv Brand. 100/PK



Additional Info

First hydrolytic, Class A, pharmaceutical grade, Borosilicate Glass.

Use with 8mm aluminum crimp caps, 8mm or 9mm push caps. These flat bottom vials have serum finish tops that provide the best seal for chromatography, when used properly with correct aluminum crimp caps.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$37.58

Printed On 2024/05/04