

Catalog Number

9509C-0CV-AQ

Vials, Snap Top, Glass. Clear, 2ml. An 11mm snap-ring and 12x32mm outer dimensions. For use as an autosampler vial. AQ Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type-33 Borosilicate Glass.

Compares to HPLC "Certified" vials for performance but lower priced and a great value. Bottom depths and heights are designed, and QC controlled for maximum recovery and needle safety.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$30.50

Printed On 2024/05/01