

Catalog Number

9509C-3WAP-AQ

Vial & Cap Kit includes 100 EA 2ml, Snap Top, Amber Autosampler Vials with Write-On Patch and fill lines and 100 EA matching Snap Caps with Tan Silicone Rubber / White PTFE Pre-Slit Septa fitted in the Blue Caps.



Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps. The kit packaging will open and re-close making it useful for keeping vial and cap types matched, organized and clean during use.

Click on links below for more Specifications for the 11mm Vials and Caps:

Supplied with Vials: Cat. No [9509C-WAV-AQ](#) and Caps: Cat No. [9502C-31-MB](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$103.30

Printed On 2024/04/30