

## Catalog Number

9509C-3WAP-RS

**Vial & Cap Kit** includes 100 EA 2ml, Snap Top, Amber RSA™ Autosampler Vials with Write-On Patch and fill lines and 100 EA matching Snap Caps with Silicone Rubber / PTFE Pre-Slit, ultra-pure



### Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps. The kit packaging will open and re-close making it useful for keeping vial and cap types matched, organized and clean during use.

Click on links below for more Specifications for the 11mm Vials and Caps:

Supplied with Reduced Surface Activity Vials: Cat. No [9509C-WAV-RS](#) and Caps: Cat No. [9509C-30CB](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$0.00**

Printed On 2024/05/04