


Catalog Number

9509C-50CB

Caps, Snap Top, with PTFE, ultra-pure, 2mil Septa, fitted in smooth, polyethylene, Blue Caps. For use with 11mm, 12x32mm, Snap Top autosampler vials. AQR Brand. 100/PK.



Additional Info

This PTFE (only) Septa is for single injections or non-volatile samples as it will not reseal.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$0.00

Printed On 2024/05/07