

Catalog Number

9509C-AQR-3S

Caps, Crimp Top, Silicone Rubber / PTFE, ultra-pure, Pre-Slit Septa, fitted in aluminum, Silver Caps. For use with 11mm Crimp or Snap Top, 12x32mm autosampler vials. AQR Brand. 100/PK.



MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$0.00

Printed On 2024/04/29