


Catalog Number

9509C-WAV-RS

Vials, Snap Top, Glass. Amber, 2ml, with Write-On Patch and fill lines. An 11mm snap-ring and 12x32mm outer dimensions. For use as an autosampler vial. RSA Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type N51A, Borosilicate Glass, amberized. Write-On Patch and fill lines are ceramic and will resist smearing or scratching when working with solvents.

Bottom depths and heights are designed, and QC controlled for maximum recovery and needle safety. Low adsorbing and will not add Na to sample solutions or compounds that cause pH changes.

Click [HERE](#) for more information about RSA autosampler vials.

Price
\$0.00

Additional Info

This is a popular 2ml amber Snap Top autosampler vial because of the consistency of the glass chemistry from batch to batch, vial to vial. Reduced surface activity is tested for each lot and certified. Excellent for low abundance analytes, pH sensitive analytes and for LCMS.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$0.00