


## Catalog Number

9509S-10B

**Caps, Screw Top**, with White Silicone Rubber / White PTFE, ultra-low bleed Septa, bonded in "knurled", polypropylene, Dark Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. AQ Brand.



### Additional Info

This very high purity Septa has low extractables and is excellent choice for HPLC. Compares to "Certified" Brands.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$35.19**

Printed On 2024/07/27