

Catalog Number

9509S-1WCP-RS-A

Vial & Cap Kit Includes 100 EA 2ml, Screw Top, Clear RSA™ Autosampler Vials with Write-On Patch and fill lines and 100 EA matching Screw Caps with Clear AQR Silicone Rubber / Clear PTFE,



Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps.

Click on links below for more Specifications for the 9mm Vials and Caps:

Supplied with Reduced Surface Activity Vials: Cat. No [9509S-WCV-RS](#) and LCMS/MS compatible Caps: Cat No. [9509S-10C-B](#).

Price
\$0.00

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$0.00