

Catalog Number

9509S-3CP-AQA

Vial & Cap Kit Includes 100 EA 2ml, Screw Top, Clear Autosampler Vials and 100 EA matching Screw Caps with Tan Silicone Rubber / White PTFE, Pre-Slit Septa bonded in the Blue Caps. Vials and Caps are Pre-



Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps.

Click on links below for more Specifications for the 9mm Vials and Caps:

Supplied with Vials: Cat. No. [9509S-0CV-AQ](#) and Caps: Cat No. [9502S-31-MB](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$95.51

Printed On 2024/05/04