

Catalog Number

9509S-WAV-AQ

Vials, Screw Top, Glass. Amber, 2ml, with Write-On Patch and fill lines, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. AQ Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type N51A, Borosilicate Glass, amberized. Write-On Patch and fill lines are ceramic and will resist smearing or scratching when working with solvents.

Compares to HPLC "Certified" vials for performance but lower priced and a great value. Bottom depths and heights are designed, and QC controlled for maximum recovery and needle safety.

Amber vials are often used for light sensitive compounds and considered low actinic.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$40.13

Printed On 2024/04/30