

Catalog Number

9509S-WCV-AQ

Vials, Screw Top, Glass. Clear, 2ml, with Write-On Patch and fill lines, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. AQ Brand. 100/PK.



Additional Info

Upgrade Your Lab with Our Premium Glass Vials | First Hydrolytic, Class A, Type-33 Borosilicate Glass | Durable Ceramic Write-On Patch | High-Quality Alternative to Waters Certified Brand number 186000273

Why Choose Our Vials?

Our vials are made from first hydrolytic, Class A, Type-33 borosilicate glass, ensuring superior quality and performance. The ceramic write-on patch and fill lines resist smearing or scratching, even when working with solvents, providing durability and

Price
\$30.96

Additional Info

ease of use. Save money and do not sacrifice manufacturers "compliance" to their warranty specifications or vial performance.

Key Benefits:

- High-Quality Material: First hydrolytic, Class A, Type-33 borosilicate glass for exceptional performance.
- Durable Write-On Patch: Ceramic write-on patch and fill lines resist smearing and scratching.

Cost-Effective and Reliable:

- Matches the performance of Waters "Certified" and other vials but at a lower price, offering great value.
- Designed for Sample and Instrument Safety: Bottom depths and heights are meticulously designed and QC controlled for maximum recovery and needle safety.

Optimize Your Lab's Efficiency:

- Maximum Recovery: Ensures you get the most out of your samples.
- Needle Safety: Designed to protect your equipment and ensure accurate results.

Upgrade Your Lab Supply Chain Today!

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$30.96