

## Catalog Number

9512C-0DCV

**Vials, Snap Top, Glass. Clear, 1.8ml, Max Recovery™. Deactivated. An 11mm snap-ring and 12x32mm outer dimensions. For use as an autosampler vial. AQ Brand. 100/PK.**



### Additional Info

First hydrolytic, Class A, Type-33 Deactivated / Silanized, Borosilicate Glass. This 1.8ml vial has the low residual volume of 8ul due to the center-draining effects of smooth walls (no ledges) and a conical bottom, which allows needles with bottom entry ports to recover the maximum amount during use.

Bottom depths are designed and QC controlled for the most critical autosampler to prevent needles hitting the bottom of the vials.

*This vial has a large headspace with conical bottom for maximum recovery of sample. Silanized to deactivate glass surface for biological and other problematic compounds.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$258.29**

Printed On 2024/05/07