

## Catalog Number

9512C-0E-0CV-H

**Vials, Crimp Top, Glass.** Clear, 1.8ml, Max Recovery Style. An 11mm Serum Finish with 12x32mm outer dimensions and rounded bottom. For use as an autosampler vial. MicroSolv Brand. 100/PK.



### Additional Info

First hydrolytic, Class A, Type N51A, Borosilicate Glass. Features an 25ul residual volume for most autosamplers due to the center draining interior and conical bottom. Rounded exterior bottom helps when inserting vials into racks or autosamplers with tight fittings.

Use with 11mm aluminum crimp caps. These lower cost, center draining vials have serum finish tops that provide the best seal for chromatography, when used properly with correct aluminum crimp caps.

*Lowest cost in class for budget minded laboratories.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$70.64**

Printed On 2024/05/02