

Catalog Number

9512C-3CP-RS

Vial & Cap Kit includes 100 EA 1.8ml, Maximum Recovery, Snap Top, Clear RSA™ Autosampler Vials and 100 EA matching Snap Caps with Clear AQR Silicone Rubber / Clear PTFE Pre-Slit, ultra-pure Septa, fitted



Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps. The kit packaging will open and re-close making it useful for keeping vial and cap types matched, organized and clean during use.

Click on links below for more Specifications for the 11mm Vials and Caps:

Supplied with Reduced Surface Activity Vials: Cat. No [9512C-0CV-RS](#) and Caps: Cat No. [9509C-30CB](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$275.98

Printed On 2024/05/05