

## Catalog Number

9512C-AQR-1CP

**Vial & Cap Kit** includes 100 EA 1.2ml, MRQ, Snap / Crimp Top, Clear Autosampler Vials and 100 EA matching Crimp Caps with Clear AQR Silicone Rubber / Clear PTFE, ultra-pure Septa, fitted in the Silver Caps.



### Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps. The kit packaging will open and re-close making it useful for keeping vial and cap types matched, organized and clean during use.

Click on links below for more Specifications for the Crimp Top Vials and Caps:

Supplied with Vials: Cat. No [9512C-0CV-T](#) and Caps: Cat No. [9509C-AQR-1S](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$223.34**

Printed On 2024/05/02