

Catalog Number

9512S-0AV

Vials, Screw Top, Glass. Amber, 1.8ml, Max Recovery™, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. AQ Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type N51A pharmaceutical grade, Borosilicate Glass, amberized. Amber glass is often used for light sensitive compounds and is considered low actinic.

This market leading 1.8ml vial has the low residual volume of 8ul due to the center-draining effects of smooth walls (no ledges) and a conical bottom, which allows needles with bottom entry ports to recover the maximum amount during use.

Bottom depths are designed and QC controlled for the most critical autosampler to

Price**\$304.46**

Additional Info

prevent needles hitting the bottom of the vials. This product is manufactured in an ISO 9001:2015 Certified facility.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$304.46